

Face Shields

Technical Data Sheet

KEY FEATURES

- Lightweight, anti-fog protection from splashing or spraying
- Foam Headband with Elastic Strap
- Can be used in conjunction with a N95 mask, goggles, glasses or be worn alone without fogging up
- Made from polyethylene terephthalate (PET)

Color	Clear Shield / Blue Headband
Dimensions	9 Inch Length Shield / 1 X 1 X 9 Inch Foam Strip
Fastening Type	Foam Headband / Elastic Strap
Lens Tint	Clear Tint
Size	One Size Fits Most
Sterility	NonSterile
Style	Full Length
Style	Anti-fog



APPROVALS AND STANDARDS

- Compliance to ANSI Z87.1 -2015 Occupational and Educational Personal Eye and Face Protection Devices
- Clause 5.1.2 - Luminous Transmission
- Clause 8.1.2 - Droplet and Splash Hazard

ORDERING INFORMATION

SKU #	DELIVERY	SHIELDS/BOX	BOXES/CASE	SHIELDS/CASE
J2FSMG01B	7-28 Days	12	20	240

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.