

# Nitrile Examination Gloves

## Technical Data Sheet

### MANUFACTURER

- Intco (United States: Basic)

### MATERIAL & DESIGN

- Nitrile Synthetic, latex free, powder free
- Beaded cuff for easy donning
- Microtextured fingertips and fully textured exterior surface for improved wet and dry grip
- Stretchable synthetic nitrile material for comfort
- Color - Blue

### THICKNESS

- 4 ml
- Fingers 0.12 mm
- Palms 0.08

### SIZE

- Small, Medium, Large, Extra-Large



### APPROVALS AND STANDARDS

- FDA Approval- 510(k) Approval Number: K990186
- Meets ASTM Standard ASTM D3578-95

## ORDERING INFORMATION

	SKU #	DELIVERY	GLOVES/BOX	BOX/CASE	GLOVES/CASE
#1	J2NG4MSMC (SM)	14-32 Days	100	10	1,000
#2	J2NG4MMDC (MD)				
#3	J2NG4MLGC (LG)				
#4	J2NG4MXLC (XL)				

### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.