


# AAMI Level 3 Disposable Isolation Gowns - ("New York")

## Technical Data Sheet

### KEY FEATURES

- Full back coverage
- Elastic on the wrist for comfortable fit
- 50gsm Polypropylene Non-woven
- Provides comfort while providing protection from moderate to high fluid exposure
- 50gsm Polypropylene Non-woven
- Sizes - One Size Fits Most
- Product Color: Blue



### AAMI PB70 Guidelines

	LEVEL 1	LEVEL 2		LEVEL 3	
Test	AATCC 42:2000	AATCC 42:2000	AATCC 127:1998	AATCC 42:2000	AATCC 127:1998
Requirements at 4% AQL	Water Impact ≤ 4.5 g	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 20 cm	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 50 cm
Anticipated Risk of Exposure	Low	Moderate		Moderate to High	

For additional information on guidelines on gowns, please refer to the following FDA link - <https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns>

#### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

# AAMI Level 3 Disposable Isolation Gowns - ("New York")

## Technical Data Sheet - (Continued)

### INTENDED USAGE

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.



### ORDERING INFORMATION

SKU #	DELIVERY	BOX/CASE COUNT
J2IG3SSMC	7-28 Days	100 units per case

#### For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.

**J2** | MEDICAL  
SUPPLY  
J2MedicalSupply.com