

AAMI Level 2 Disposable Isolation Gowns - ("Ocean")

Technical Data Sheet

KEY FEATURES

- Single layer polyethylene
- Reliable protection from transfer of microorganisms and body fluids
- Over the head with thumb loop
- Perforated back for easy removal
- Waist and neck tie closures
- Lightweight, breathable and comfortable
- Latex Free
- Sizes - One Size Fits Most
- Product Color: Blue



AAMI PB70 Guidelines

| | LEVEL 1 | LEVEL 2 | | LEVEL 3 | |
|------------------------------|----------------------|----------------------|------------------------------|----------------------|------------------------------|
| Test | AATCC 42:2000 | AATCC 42:2000 | AATCC 127:1998 | AATCC 42:2000 | AATCC 127:1998 |
| Requirements at 4% AQL | Water Impact ≤ 4.5 g | Spray Impact ≤ 1.0 g | Hydrostatic Pressure ≥ 20 cm | Spray Impact ≤ 1.0 g | Hydrostatic Pressure ≥ 50 cm |
| Anticipated Risk of Exposure | Low | Moderate | | Moderate to High | |

For additional information on guidelines on gowns, please refer to the following FDA link - <https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns>

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

AAMI Level 2 Disposable Isolation Gowns - ("Ocean")

Technical Data Sheet - (Continued)

INTENDED USAGE

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.



ORDERING INFORMATION

| SKU # | DELIVERY | BOX/CASE COUNT |
|-----------|-----------|--------------------|
| J2IG20CNC | 7-28 Days | 100 units per case |

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.

J2 | MEDICAL
SUPPLY
J2MedicalSupply.com