


AAMI Level 2 Disposable Isolation Gowns - ("Valencia")

Technical Data Sheet

KEY FEATURES

- Over the head with thumb loop
- Functionally designed for easy donning and removal
- Available in four sizes - SM/MD, LG/XL, 2X, and 3X
- 70g Spunlace nonwoven
- Superior comfort while providing protection from light to moderate fluid exposure
- Product Color: Blue



AAMI PB70 Guidelines

	LEVEL 1	LEVEL 2		LEVEL 3	
Test	AATCC 42:2000	AATCC 42:2000	AATCC 127:1998	AATCC 42:2000	AATCC 127:1998
Requirements at 4% AQL	Water Impact ≤ 4.5 g	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 20 cm	Spray Impact ≤ 1.0 g	Hydrostatic Pressure ≥ 50 cm
Anticipated Risk of Exposure	Low	Moderate		Moderate to High	

For additional information on guidelines on gowns, please refer to the following FDA link - <https://www.fda.gov/medical-devices/personal-protective-equipment-infection-control/medical-gowns>

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

AAMI Level 2 Disposable Isolation Gowns - ("Valencia")

Technical Data Sheet - (Continued)

INTENDED USAGE

- Single-use gown to protect the wearer from the transfer of microorganisms and body fluids in low or minimal risk patient isolation situations.
- For use in patient care areas and other points of interaction where there is a need for barrier protection.
- Not intended for high risk surgical procedures, or where the risk of contamination is high.



ORDERING INFORMATION

SKU #	DELIVERY	BOX/CASE COUNT
#1 J2IG2VSMC (SM/MD)	7-28 Days	80 units per case
#2 J2IG2VLXC (LG/XL)		
#3 J2IG2V2XC (2X)		
#4 J2IG2V3XC (3X)		

For More Information:

855.615.8633

Orders@J2MedicalSupply.com

© J2 Medical Supply 2020. All rights reserved.

J2 | MEDICAL
SUPPLY
J2MedicalSupply.com